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Uso Diagnostico in Vitro / In Vitro Diagnostic use only

PRINCIPLE

HbA1c Calibration Set is a set of 5 calibrators consisting of human erythrocytes. It can be used for the quantitative determination of human glycosylated hemoglobin (A1c) with turbidimetric method.

Exact values are indicated in the value sheet.

REAGENT COMPOSITION

Calibrators L1-L5: Human erythrocytes. - Stabilizers

REAGENT PREPARATION, STORAGE AND STABILITY

Calibrators L1-L5 are lyophilized. Reconstitute with 1 mL of H₂O. Shake gently avoiding foaming. Wait 10 minutes before use.

Calibrators are ready to use after reconstitution. Hemolyzing step is not required. Stability is 30 days at 2-8° C. Do not use the reagents after the expiration date.

PRECAUTIONS AND WARNINGS

The products do not contain dangerous substances or mixtures, according to the EC Regulation n° 1272/2008 or their concentrations are such as not to be considered persistent, bioaccumulative or toxic (PBT). The product is classified and labeled in accordance with EC directives or respective national laws. Sodium azide, less than 0.1%. However, in compliance with the normal prudential rules that everyone has to keep managing any chemical or laboratory reagent, in case of contact of reagents with the operator, you must apply the following first aid:

S26 (P305 – P351 – P338): In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 (P302 – P352): After contact with skin, wash immediately with plenty of water.

S36/37/39 (P280): Wear suitable protective clothing, gloves and eye/face protection.

S46 (P301 – P310): If swallowed, seek medical advice immediately and show container or label. If victim is conscious and alert, wash out mouth with water.

S56 (P273): Dispose of this material and its container at hazardous or special waste collection point.

S63 (P304 – P340): In case of accident by inhalation: remove casualty to fresh air and keep at rest. If breathing is difficult, give oxygen and get medical aid.

The human components present in the mixture were tested with CE-approved and FDA-approved methods and found to be negative for HBsAg, HIV 1 & 2, HCV.

However, since there are no methods to ensure the complete absence of such infectious agents, this product must be handled as if it were risky and potentially capable of transmitting infectious diseases of any kind, in accordance with Good Laboratory Practice standards.

FOR MORE INFORMATION, REQUEST SAFETY DATA SHEET FOR REAGENT (MSDS) AT THE MANUFACTURER.

SIMBOLOGY

CE	CE Mark (EC Directive 98/79)		
	In Vitro Diagnostic		Temperature Limitation
	Consult instructions for use		Contains sufficient for <n> test
	Catalog Number		Use By
	Batch Code		Manufacturer

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