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### PRINCIPLE

The D-Dimer Control Set Master is a set of two serum levels. These two levels are to be used exclusively with the D-Dimer Assay which is a latex immunoturbidimetric assay for the quantification of D-Dimer in human plasma.

### REAGENT COMPOSITION

#### Reactive components

- Human serum and additives

#### Non-reactive components

- Sodium azide (NaN<sub>3</sub>) <0.1%

### REAGENT PREPARATION, STORAGE AND STABILITY

The D-Dimer Control Set Master is a two-level set that is supplied in lyophilized form (2x1 mL). The controls are manufactured from human serum. The concentrations of the controls are lot specific and D-Dimer concentrations are expressed in µg/mL FEU. The D-Dimer Control Set Master is supplied in lyophilized form. Add 1.0 mL of distilled water to the vial. Cap and leave at room temperature for 30 minutes, inverting gently to dissolve contents. Mix carefully before each use. Store at 2 to 8 °C. DO NOT FREEZE.

Note: Store controls tightly capped when not used.

#### Stability

- Unopened vials: stability is up to the expiration date on label when stored at 2 to 8 °C
- After reconstitution: stability is 14 days when store at 2 to 8 °C

### PRECAUTIONS AND WARNINGS

The products do not contain dangerous substances or mixtures, according to the EC Regulation n° 1272/2008 or their concentrations are such as not to be considered persistent, bioaccumulative or toxic (PBT). The product is classified and labeled in accordance with EC directives or respective national laws. Sodium azide, less than 0.1%. However, in compliance with the normal prudential rules that everyone has to keep managing any chemical or laboratory reagent, in case of contact of reagents with the operator, you must apply the following first aid:

S26 (P305 – P351 – P338): In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 (P302 – P352): After contact with skin, wash immediately with plenty of water.

S36/37/39 (P280): Wear suitable protective clothing, gloves and eye/face protection.

S46 (P301 – P310): If swallowed, seek medical advice immediately and show container or label. If victim is conscious and alert, wash out mouth with water.

S56 (P273): Dispose of this ial and its container at hazardous or special waste collection point.

S63 (P304 – P340): In case of accident by inhalation: remove casualty to fresh air and keep at rest. If breathing is difficult, give oxygen and get medical aid.

The human components present in the mixture were tested with CE-approved and FDA-approved methods and found to be negative for HBsAg, HIV 1 & 2, HCV.

However, since there are no methods to ensure the complete absence of such infectious agents, this product must be handled as if it were risky and potentially capable of transmitting infectious diseases of any kind, in accordance with Good Laboratory Practice standards.

FOR MORE INFORMATION, REQUEST SAFETY DATA SHEET FOR REAGENT (MSDS) AT THE MANUFACTURER.

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	CE Mark (EC Directive 98/79)		
	In Vitro Diagnostic		Temperature Limitation
	Consult instructions for use		Contains sufficient for <math>N</math>-test
<b>REF</b>	Catalog Number		Use By
<b>LOT</b>	Batch Code		Manufacturer

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