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ASO CALIBRATOR SET is a set of liquid calibration sera joined to biological additives and antibacterial agents and then diluted in saline solution. The constituent concentrations and activities have been selected to ensure an optimum control for use in manual and automatic analyzers. Components' concentrations are lot-specific.

Exact values are indicated in the value sheet.

PREPARATION AND STABILITY

R1 to R4: Liquid and ready to use.

Shortly before use, swirl gently several times avoiding formation of foam (inverting instead of shaking the vial).

Stability:

Unopened vial: up to expiration date at 2-8 °C
 Once vial is opened: 1 month at 2-8 °C
 Onboard stability: 24 hours

PRECAUTIONS AND WARNINGS

The products do not contain dangerous substances or mixtures, according to the EC Regulation n° 1272/2008 or their concentrations are such as not to be considered persistent, bioaccumulative or toxic (PBT). The product is classified and labeled in accordance with EC directives or respective national laws. Sodium azide, less than 0.1%.

However, in compliance with the normal prudential rules that everyone has to keep managing any chemical or laboratory reagent, in case of contact of reagents with the operator, you must apply the following first aid:

S26 (P305 – P351 – P338): In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 (P302 – P352): After contact with skin, wash immediately with plenty of water.

S36/37/39 (P280): Wear suitable protective clothing, gloves and eye/face protection.

S46 (P301 – P310): If swallowed, seek medical advice immediately and show container or label. If victim is conscious and alert, wash out mouth with water.

S56 (P273): Dispose of this material and its container at hazardous or special waste collection point.

S63 (P304 – P340): In case of accident by inhalation: remove casualty to fresh air and keep at rest. If breathing is difficult, give oxygen and get medical aid.

The human sera present in the mixture were tested with CE-approved and FDA-approved methods and found to be negative for HBsAg, HIV 1 & 2, HCV. However, since there are no methods to ensure the complete absence of such infectious agents, this product must be handled as if it were risky and potentially capable of transmitting infectious diseases of any kind, in accordance with Good Laboratory Practice standards.

FOR MORE INFORMATION, REQUEST SAFETY DATA SHEET FOR REAGENT (MSDS) AT THE MANUFACTURER.

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SIMBOLOGY

	CE Mark (EC Directive 98/79)		
	In Vitro Diagnostic		Temperature Limitation
	Consult instructions for use		Contains sufficient for <n> test
REF	Catalog Number		Use By
LOT	Batch Code		Manufacturer

only

TABELLA VALORI / VALUE TABLE (Unità U/L / Units IU/L)

R1

R2

R3

R4

* Titolato contro Standard WHO / Titled against WHO Standardization

** Titolato contro Standard Siemens / Titled against Siemens Standardization