

Dia-CaCl₂



CALCIUM CHLORIDE SOLUTION FOR COAGULATION TESTS

Cat. No.: 41192 12 x 16 ml
 Cat. No.: 41048 12 x 4 ml

- Different levels of control for quality control (Dia-CONT I-II; Cat. No.: 91020, 91010).
- Optical or mechanical coagulation analyser for measuring, Diagon analysers (Coag Line) are recommended.

PRODUCT NAME
 Dia-CaCl₂ solution.

INTENDED USE
(For In Vitro Diagnostic Use Only)
 Dia-CaCl₂ is a supplementary reagent of a laboratory, two steps screening test of APTT.

DIAGON LTD.
 Baross u. 48-52, 1047 Budapest, Hungary
 Tel.: +36 1 3696500
 Fax.: +36 1 3696301
 Web: www.diagon.com
 E-mail: diagon@diagon.com

ACTIVE INGREDIENTS
 Dia-CaCl₂ is a 0,025 M buffered solution with preservative.

- PRECAUTIONS**
- Person installing the Dia-CaCl₂ solution must be a trained laboratory professional!
 - Reagent coming into contact with specimens and other materials should be handled as if capable of transmitting infection and should be disposed of with proper precautions!
 - Avoid microbial contamination of the reagent otherwise erroneous results may occur!
 - All reagents, waste and utilized disposable laboratory equipment should be considered as hazardous waste! Their handling and disposal should be done according to the valid hazardous material processing regulation.
 - Do not use the reagent beyond the expiration date printed on the label!

SYMBOLS			
	Manufacturer		Use-by date
	Batch code		Catalogue number
	Do not use if package is damaged		Fragile, handle with care
	Keep dry		Temperature limit
	Caution		Consult instruction for use
	Contains sufficient for <n> tests		In vitro diagnostic medical device
	CE mark		This side up

PREPARATION
 Dia-CaCl₂ solution is ready for use. Wait until the reagent reaches the working temperature!

TEST PROCEDURE
 Follow as described in the instruction for use of Dia-PTT or Dia-PTT LIQUID tests to measure the required parameter.

STORAGE AND STABILITY
 Dia-CaCl₂ solution in intact vial is stable until the expiration date given on the vial, when stored at 2-8°C. Stability after opening in the original vial is shown in below table:

T (°C)	37	15-19	2-8
Week	1 day	8	8

- MATERIALS REQUIRED BUT NOT PROVIDED**
- APTT:
 - Dia-PTT (Cat. No.: 71024; 71030; 71048)
 - Dia-PTT LIQUID (Cat. No.: 72012; 72024; 72048; 72096)