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BIBLIOGRAFIA

- WHO Expert Committee on Standardisation: WHO Technical Report Series. No.889; 1999.

CONTROL PLASMAS FOR COAGULATION TESTING – 2 LEVELS

REF CO1110 7+7x1 mL

In Vitro Diagnostic Use

INTENDED USE

PLASMA CONT I-II is a set of lyophilized control plasmas, at 2 levels (normal and abnormal) for the Internal Quality Control (Precision) of the main coagulation tests (PT, PTT, Fibrinogen, Antithrombin). To verify Accuracy, refer to EQC (External Quality Control) programs managed by Certified Authority.

The application of an Internal Quality Control system is strongly recommended to establish the repeatability of the analytical data.

REAGENT COMPOSITION

PLASMA CONT I-II: Pool of lyophilized human plasmas.

REAGENT PREPARATION AND STABILITY

PLASMA CONT I-II: Lyophilized. Before reconstitution, if stored at 2-8 °C, it is stable until the expiration date shown on the label. Reconstitute with 1 mL of distilled water. Keep the reconstituted reagent 30 minutes at 18-25° C before using it, gently shaking it, avoiding the liquid touching the cap and foaming.

Stability after reconstitution: 4 hours at 22°C, 30 days at ≤ -20°C.

Freeze and thaw only once.

PROCEDURE

After reconstitution of the plasmas, use them as samples under examination. Record the results obtained and calculate the mathematical control parameters (Average, SD, CV%) to check the Precision. Using charts, allows a better evaluation of daily data.

Compare the results obtained with the specific table for the lot enclosed in the kit.

LIMITATIONS

- Use the attached table of values as a guideline. The central value of the reference range does not necessarily have to be obtained.
- The product lends itself to being used as a verification of the level of Precision (Repeatability).
- To check the level of Accuracy (Reproducibility), use a EQC service (External Quality Control) managed by Certified Authority.
- Each laboratory should determine its own reference range referring its instrumentation and handling.

PRECAUTIONS

The human sera present in the mixture were tested with CE-marked and FDA-approved methods and were found to be negative for the presence of HBsAg, HIV 1 & 2, HCV. However, since there are no methods capable of ensuring the total absence of such infectious agents, this product must be handled as if it were risky and potentially capable of transmitting infectious diseases of any kind, in accordance with the rules of Good Laboratory Practice.

The product does not contain dangerous substances, according to the EC Regulation n ° 1272/2008 or their concentrations are such as not to be considered persistent, bioaccumulative or toxic (PBT). The product is classified and labeled in accordance with EC directives or respective national laws. Sodium azide, less than 0.1%.

However, in compliance with the normal prudential rules that everyone has to keep managing any chemical or laboratory reagent, in case of contact of reagents with the operator, you must apply the following first aid:

S26 (P305 – P351 – P338): In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 (P302 – P352): After contact with skin, wash immediately with plenty of water.

S36/37/39 (P280): Wear suitable protective clothing, gloves and eye/face protection.

S46 (P301 – P310): If swallowed, seek medical advice immediately and show container or label. If victim is conscious and alert, wash out mouth with water.

S56 (P273): Dispose of this material and its container at hazardous or special waste collection point.

S63 (P304 – P340): In case of accident by inhalation: remove casualty to fresh air and keep at rest. If breathing is difficult, give oxygen and get medical aid.

All test specimens, calibrators and controls should be treated as potentially infectious material capable of transmitting HIV, hepatitis or any other viral infection not currently known.

FOR MORE INFORMATION, REQUEST SAFETY DATA SHEET FOR REAGENT (MSDS) AT THE MANUFACTURER.

SIMBOLOGY

CE	CE Mark (EC Directive 98/79)		
	In Vitro Diagnostic		Temperature Limitation
	Consult instructions for use		Contains sufficient for <n> test
	Catalog Number		Use By
	Batch Code		Manufacturer

BIBLIOGRAFY

- 1 WHO Expert Committee on Standardisation: WHO Technical Report Series. No.889; 1999.