PRODUCT IDENTIFICATION:

Name: DILUENT; SHEATH; LYSE; RINSE;
 PROBE CLEANER 3 DIFF H3;
 PROBE CLEANSER 5 DIFF H5

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INTENDED USE

For In Vitro Diagnostic use only. Reagent to be used with cell counters

CONSERVATION AND STABILITY AFTER FIRST OPENING

Keep between +15/+30°C. Install every reagent only after verifying integrity of the kit and validity of device (expiry date). Once reagent is open, it can be used for consecutive 60 days. In case of temporary suspension, use the intact original cap to close the container.

COMPOSITION

Reagents		Active ingredients	
Product Name	REF	Name	Content
DILUENT	M3001	sodium chloride	< 5 g/l
	M5001	sodium sulfate	< 12 g/l
SHEATH	M5003	aromatic alcohol	< 1 g/l
LYSE	M3003	quaternary ammonium	< 5 g/l
	M5004	salts	
RINSE	M3005	surfactant	< 3 g/l
	M3002		
PROBE CLEANER 3 DIFF H3	M3004	hypochlorite	< 3 %
PROBE CLEANSER 5 DIFF H5	M5005	hypochlorite	< 3 %

WARNINGS

Reagents		Info on Labelling
Product Name	REF	
DILUENT	M3001	Danger label under Regulation (EC) 1272/2008 CLP not required
	M5001	
SHEATH	M5003	Danger label under Regulation (EC) 1272/2008 CLP not required
LYSE	M3003	Danger label under Regulation (EC) 1272/2008 CLP not required
	M5004	
RINSE	M3005	Danger label under Regulation (EC) 1272/2008 CLP not required
	M3002	
PROBE CLEANER 3 DIFF H3	M3004	Danger label under Regulation (EC) 1272/2008 CLP not required
PROBE CLEANSER 5 DIFF H5	M5005	Danger label under Regulation (EC) 1272/2008 CLP not required

DISPOSAL INFORMATION

The remains of the products and waste from their use as well as containers must be disposed of in accordance with current legislation. Waste of normal use and potential infected substances must be carefully manipulated and disposed of in accordance with current legislation as well.

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PROCEDURE

DILUENT is a ready to use isotonic solution (diluent).
SHEATH is a ready to use lyser solution favoring WBC transport.
LYSE is a ready to use lyser solution for hemoglobin.
RINSE is a ready to use detergent solution/zeroing hemoglobin.
PROBE CLEANER 3 DIFF H3 is a cleaning solution for needles and tubes for 3 diff cell counters.
PROBE CLEANER 5 DIFF H5 H3 is a cleaning solution for needles and tubes of 3 and 5 diff cell counters.

INSTALLATION MODE : for the installation of reagents, the start-up of the lines and how to operate the instrument refer to the Operating Manual of the instrument in use. After each installation of a new batch, check the bottom (background) following the User Manual of the instrument. Operate at the temperature indicated in the User Manual. Do not transfer remnants from one container to another : operate according to good laboratory practice (GLP) in order to avoid any kind of contamination both of the external environment towards the device and from the device to the external environment.

When installing reagents for the first time, it is recommended to empty the lines and rinse with deionized water in order to avoid any kind of contamination of the reagent that could compromise the performance of the analytical system (instrument). After the initial start-up of the line, background counts must be made to verify that results are within the limits specified by the Operating Manual. Check on instrument calibration.

ANALYSIS METHOD AND CALCULATION OF RESULTS: the information on the principles of analysis on which the determination of the different analytical parameters is based and how the results are calculated are to be considered on consultation of the User Manual of the instrument in use.

INTERNAL QUALITY CONTROL

- Background verification : After filling instrument with the new reagents, perform a background check count following the instructions as reported in the Instrument's User Manual. Make sure that the resulting values, for each measured parameter, fall within the limits specified in the Manual. If one or more parameters do not fall within the specified values, check for absence of bubbles in the circuits, and follow instructions in accordance with the requirements of the manufacturer of the instrument. If the problem is attributable to one of the reagents, replace the package with an integral one and repeat the verification after filling the load line/circuit of the reagent. Do not proceed with the analyses until the values are acceptable.

- Calibration check : it is recommended to check the calibration of the instrument at each replacement even of a single reagent and in any case after each ordinary and extraordinary maintenance and in case of a significant change (drift) of accuracy or if results do not fall within the limits of acceptability predetermined in accordance with the instructions in the User Manual of the instrument.

Before performing calibration procedures, it is recommended that you perform preventive maintenance and verify accuracy and precision, as per the User Manual. Use appropriate material such as fresh whole blood whose values have been assigned using manual reference methods or with appropriate calibrator following the instructions of the instrument's User Manual.

- Daily quality check : to ensure the accuracy and reliability of the results obtained, it is recommended to perform the quality check in accordance with the procedures indicated in the User Manual of the instrument using appropriate control materials (fresh blood or stabilized hematological control) : use at least two levels (abnormal low and normal) to be analyzed at least twice before each analytical sequence and repeat the Normal level even during and at the end of each analytical sequence and after the end of a new analytical sequence on change of reagent LOT.