CONTROLLI PER VITAMINA D SIERICA E PLASMATICA/ CONTROLS FOR SERIC AND PLASMATIC VITAMIN D

REF TUC1090 2x1 mL

LOT XXXXXXXX

DD.MM.YYYY

Uso Diagnostico in Vitro / In Vitro Diagnostic use only

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VITAMIN D CONTROL SET 2 LEVELS is a set of liquid controls (2 levels) obtained from defibrinated human sera diluted in saline solution. The concentration contained in each bottle have been selected so as to ensure an optimal control both in manual use and with automatic systems. Concentrations are lot specific. **Expected or target values are indicated in the value sheet.**

PREPARATION AND STABILITY

R1 / R2: Ready to use.

Shortly before use, swirl gently several times avoiding formation of foam (inverting instead of shaking the vial).

Stability:

Unopened vial: up to expiration date at 2-8 °C

Once vial is opened: 1 month at 2-8 °C

Onboard stability: 24 hours

PRECAUTIONS AND WARNINGS

The products do not contain dangerous substances or mixtures, according to the EC Regulation n ° 1272/2008 or their concentrations are such as not to be considered persistent, bioaccumulative or toxic (PBT). The product is classified and labeled in accordance with EC directives or respective national laws. Sodium azide, less than 0.1%.

However, in compliance with the normal prudential rules that everyone has to keep managing any chemical or laboratory reagent, in case of contact of reagents with the operator, you must apply the following first aid:

S26 (P305 – P351 – P338): In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S28 (P302 – P352): After contact with skin, wash immediately with plenty of water.

S36/37/39 (P280): Wear suitable protective clothing, gloves and eye/face protection.

S46 (P301 – P310): If swallowed, seek medical advice immediately and show container or label. If victim is conscious and alert, wash out mouth with water.

S56 (P273): Dispose of this material and its container at hazardous or special waste collection point.

S63 (P304 – P340): In case of accident by inhalation: remove casualty to fresh air and keep at rest. If breathing is difficult, give oxygen and get medical aid.

The human sera present in the mixture were tested with CE-approved and FDA-approved methods and found to be negative for HBsAg, HIV 1 & 2, HCV.

However, since there are no methods to ensure the complete absence of such infectious agents, this product must be handled as if it were risky and potentially capable of transmitting infectious diseases of any kind, in accordance with Good Laboratory Practice standards.

FOR MORE INFORMATION, REQUEST SAFETY DATA SHEET FOR REAGENT (MSDS) AT THE MANUFACTURER.

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SIMBOLOGIA

CE	CE Mark (EC Directive 98/79)		
	In Vitro Diagnostic		Temperature Limitation
	Consult instructions for use		Contains sufficient for <n> test
	Catalog Number		Use By
	Batch Code		Manufacturer

Valori di Vitamina D / Vitamin D values

Livello 1 / Level 1: yy (aaa – bbb) ng/mL

Livello 2 / Level 2: xx (aaa – bbb) ng/mL